





























































Stabilis



|  |  |  |  |  |  |  |  | |
|---|---|---|---|--|---|---|---|------|
| Cefazolin sodium | 0,5 mg/ml | PVC |  | 25°C |  | 8 |  | 2372 |
| Cefazolin sodium | 0,5 mg/ml | PVC |  | 37°C |  | 24 |  | 2372 |
| Cefazolin sodium | 0,5 mg/ml | PVC |  | 4°C |  | 14 |  | 2372 |
| Cefotaxime sodium | 0,80 mg/ml | PVC |  | 37°C |  | 6 |  | 2374 |
| Cefotaxime sodium | 0,87 mg/ml | PVC |  | 25°C |  | 24 |  | 2374 |
| Ciprofloxacin lactate | 0,025 mg/ml | PVC |  | 25°C |  | 7 |  | 4263 |
| Ciprofloxacin lactate | 0,025 mg/ml | PVC |  | 37°C |  | 48 |  | 4263 |
| Deferoxamine mesylate | 0,35 mg/ml | PVC |  | 25°C |  | 7 |  | 2156 |
| Deferoxamine mesylate | 0,35 mg/ml | PVC |  | 37°C |  | 48 |  | 2156 |
| Deferoxamine mesylate | 0,35 mg/ml | PVC |  | 4°C |  | 15 |  | 2156 |
| Ofloxacin | 0,025 mg/ml | PVC |  | 25°C |  | 7 |  | 2155 |
| Ofloxacin | 0,025 mg/ml | PVC |  | 37°C |  | 48 |  | 2155 |
| Ofloxacin | 0,025 mg/ml | PVC |  | 4°C |  | 14 |  | 2155 |



Dictionnaire

| | |
|--|---|
|  Listes récapitulatives |  Stabilité des solutions |
|  Solvant |  Dianéal® PD1 4.25% glucose (Baxter) |
|  Molécule |  Concentration |
|  Contenant |  Température |
|  Conservation |  Durée de stabilité |
|  Bibliographie |  Polyvinyl chlorure |
|  Non précisée |  Jour |
|  Heure |  A l'abri de la lumière |
|  Dictionnaire | |