

# Stabilis



							
Linezolid	0,15 -> 0,6 mg/ml			25°C		7 	2163
Linezolid	0,15 -> 0,6 mg/ml			37°C		7 	2163
Linezolid	0,15 -> 0,6 mg/ml			4°C		7 	2163
Piperacillin sodium / tazobactam	0.2 / 0.025 mg/ml			23°C		7 	2451
Piperacillin sodium / tazobactam	0.2 / 0.025 mg/ml			37°C		24 	2451
Piperacillin sodium / tazobactam	0.2 / 0.025 mg/ml			4°C		14 	2451



# Dictionnaire

 Listes récapitulatives	 Stabilité des solutions
 Solvant	 Dianéal® PD2 4.25% glucose (Baxter)
 Molécule	 Concentration
 Contenant	 Température
 Conservation	 Durée de stabilité
 Bibliographie	 Polyvinyl chlorure
 A l'abri de la lumière	 Jour
 Non précisé	 Non précisée
 Heure	 Dictionnaire