



























































# Stabilis



								
Ciprofloxacin lactate	0,025 mg/ml			20°C		42		2368
Ciprofloxacin lactate	0,025 mg/ml			37°C		42		2368
Ciprofloxacin lactate	0,025 mg/ml			4°C		42		2368
Moxifloxacin	0.025 mg/ml			25°C		7		2145
Moxifloxacin	0.025 mg/ml			37°C		3		2145
Moxifloxacin	0.025 mg/ml			4°C		14		2145
Vancomycin hydrochloride	0,025 mg/ml			20°C		28		2370
Vancomycin hydrochloride	0,025 mg/ml			37°C		7		2370
Vancomycin hydrochloride	0,025 mg/ml			4°C		28		2370



## Dictionnaire

 Listes récapitulatives	 Stabilité des solutions
 Solvant	 Dianéal® PD1 1.36% glucose (Baxter)
 Molécule	 Concentration
 Contenant	 Température
 Conservation	 Durée de stabilité
 Bibliographie	 Polyvinyl chlorure
 A l'abri de la lumière	 Jour
 Non précisée	 Dictionnaire