









































































Stabilis



							
Buprenorphine hydrochloride	0,03 mg/ml			1-9°C		180 	4314
Buprenorphine hydrochloride	0,03 mg/ml			18-24°C		180 	4314
Cefamandole nafate	285 mg/ml			25°C		72 	612
Cefamandole nafate	285 mg/ml			5°C		10 	612
Dexamethasone sodium phosphate	1 >> 4 mg/ml			22°C		28 	546
Dexamethasone sodium phosphate	1 >> 4 mg/ml			4°C		28 	546
Epoetin alfa	5000 UI/ml			30°C		84 	304
Epoetin alfa	5000 UI/ml			5°C		84 	304
Morphine sulfate	1 mg/ml			22°C		91 	115
Morphine sulfate	1 mg/ml			4°C		91 	115
Moxalactam	333 mg/ml			25°C		24 	626
Moxalactam	333 mg/ml			5°C		96 	626
Parecoxib sodium	16 & 20 mg/ml			25°C		48 	1824
Parecoxib sodium	16 & 20 mg/ml			5°C		48 	1824
Parecoxib sodium	13,3 & 20 mg/ml			25°C		48 	1824
Parecoxib sodium	13,3 & 20 mg/ml			5°C		48 	1824



Dictionnaire

 Listes récapitulatives	 Stabilité des solutions
 Solvant	 NaCl 0,9% + bactériostatique
 Molécule	 Concentration
 Contenant	 Température
 Conservation	 Durée de stabilité
 Bibliographie	 Verre
 A l'abri de la lumière	 Jour
 Non précisée	 Heure
 Non précisé	 Lumière
 Seringue polypropylène	 Dictionnaire