

Stabilis



Cafeine



Stabilité des préparations

		100 mg ® = ?	TrichoFoam® 100 g	15-30°C		180			4814
		2000 mg ® = ?	TrichoFoam® 100 g	15-30°C		180			4814



















Bibliographie

	Type	Source
4814	Revue	Polonini H, Lopes de Sousa P, Silva C, Marianni B. Compatibility of Caffeine, Clobetasol Propionate, Dutasteride, Nicotinamide, and Progesterone in TrichoFoam, a Natural Vehicle for Hair Forms. Int J Pharm Compound 2024 ;28,2:161-168



Dictionnaire

 Divers	 Forme dermique
 Stabilité des préparations	 Contenant
 Origine	 Excipient
 Température	 Conservation
 Durée de stabilité	 Biosimilaire
 Données conflictuelles	 Bibliographie
 Flacon plastique	 Poudre
 A l'abri de la lumière	 Jour
 Bibliographie	 Dictionnaire